

Special Thanks to
Dr. Hans Wilderjans


Equine **ARTHROSCOPY**

Curved Hook Knife

Curved Hook Knife

Design by Dr. Hans Wilderjans
Silicone Handle

163-200-20



Curved Hook Knife

INDICATIONS

Chronic tenosynovitis of the digital flexor tendon sheath with thickening and constriction of the palmar annular ligament.

GOAL

Section of a thickened and constricted annular ligament using a tenoscopic approach.

This knife is specially designed to facilitate easy cutting of the palmar annular ligament in cases of chronic tenosynovitis of the digital flexor tendon sheath in the horse.



SPECIFICATIONS



Hook Probe

Curette

Pick Probe

Banana Knife

Periosteal Elevator




Rita Leibinger Propylux Grip

for secure and fatigue-free working
Sterilizable up to 134°C / 273°F

Hook Probe

Propylux handle blue, 225 mm, working end 110 mm, D 18

163-600-10



Banana Knife

Curved,
Propylux handle blue, 225 mm, working end 110 mm, D 18

163-600-11



Curette

Straight
Propylux handle blue, 225 mm, working end 110 mm, D 18

Product Code	Description
163-600-12	2 mm
163-600-13	3 mm
163-600-14	4 mm
163-600-15	5 mm



Curette

Curved, 45°, 5 mm
Propylux handle blue, 225 mm, working end 110 mm, D 18

163-600-16



Pick Probe

Propylux handle blue, 225 mm, working end 110 mm, D 18

Product Code	Description
163-600-17	strongly curved
163-600-18	slightly curved



Elevator

Small
Propylux handle blue, 225 mm, D 18

Product Code	Description
163-600-19	110 mm working end
163-600-22	65 mm working end



Elevator

Large
Propylux handle blue, 300 mm, working end 150 mm, D 18

163-600-20



Periosteal Elevator

Propylux handle blue, 225 mm, working end 110 mm, D 18

163-600-21



Hook Grasping Forceps

2x2 teeth, long

163-500-28



Hook Grasping Forceps

Suprapatellar

163-501-28



Arthroscope

4.0 mm, autoclave 132°/273°F,
Stainless Steel



Product Code	Description
162-0019-00	0 °, 175 mm long
162-0019-30	30 °, 175 mm long

Arthroscopy Shaft, OD 6.0 mm

With 2 stop-cocks, for 4.0mm Arthroscopes 0°/30°

162-0040-30

Arthroscopy Obturator


For arthroscopy shaft 162-0040-30, blunt

162-0040-35





arthroscopy.leibinger.vet

Arthroscopy Rongeurs




Product Code	Description
163-404-791	Ferris-Smith Arthr. Rongeur, pike shaped tip, straight, 3x10mm, 18cm
163-404-792	Ferris-Smith Arthr. Rongeur, pike shaped tip, straight, 4x10mm, 18cm
163-404-794	Ferris-Smith Arthr. Rongeur, pike shaped tip, straight, 6x10mm, 18cm
163-404-795	Ferris-Smith Arthr. Rongeur, pike shaped tip, straight, 6x12mm, 18cm
163-404-802	Ferris-Smith Arthr. Rongeur, pike shaped tip, up-angle, 4x10mm, 18cm
163-404-812	Ferris-Smi. Arthr. Rongeur, pike shaped tip, down-angle, 4x10mm, 18cm



Pike-Shaped tip

Arthroscopy Rongeurs



Product Code	Description
163-404-755	Schlesinger Rongeur, straight, 3x10mm serrated jaw, 5" Shaft
163-404-756	Schlesinger Rongeur, straight, 4x10mm serrated jaw, 5" Shaft
163-404-768	Wilde Ethmoid Rongeur, straight, 4mm jaw, 5,5" Shaft
163-404-782	Mcllwraith Peapod Rongeur, straight, 2.5x6mm, 5,5" Shaft
163-404-784	Mcllwraith Peapod Rongeur, up-angle, medium jaw, 5,5" Shaft



Equine Arthroscopy Set

Set contains:

- Hook Probe
- Banana Knife
- Curettes 2,3,4 & 5 mm
- Curette curved, 5mm
- 2 pick probes, strongly and slightly curved
- 2 elevators, standard 110mm, large 150mm, short 65mm
- Dissector Periosteal Elevator
- 6 Ferris Smith Arthroscopic Rongeurs,
(3x10, 4x10, 6x10, 6x12, 4x10 up- & 4x10 down angle)
- 5 Arthroscopic Rongeurs,
(Schlesinger, Wilde Ethmoid, Mcllwraith Peapod)
- Hook Grasping Forceps 2x2 teeth
- Hook Grasping Forceps Suprapatellar



163-1000-00

